

# Disposable collecting cartridges

The disposable cartridges are designed for use with one- and two-liter reusable collecting containers produced by OGARIT.

The use of disposable cartridges eliminates the need to open the lid of the reusable collecting container filled with suctioned content for emptying. This completely reduces the contact of the medical staff with the container's contents which ensures the highest standards of hygiene in the workplace.



The disposable cartridges have a built-in antibacterial-hydrophobic filter which prevents the leakage of the suctioned content on the outside of the cartridge when it is fully filled-up. This allows to protect medical personnel and the environment from pathogens. For the additional protection and the convenience of use the disposable cartridges may contain a special agent, which changes the consistency of stored liquid into gel. The risk of fluid leakage out of the disposable collecting cartridge is fully reduced.

There is an additional slot in the cover allowing the easy collection of samples from the inside of the disposable cartridge.



The disposable cartridges of 1 or 2 litres capacity manufactured by OGARIT may be used in all electric aspirators and aspiration sets adapted for the central vacuum system.

## UTILIZATION

The OGARIT cartridges are disposable products. Each medical facility should provide an appropriate disposal of used cartridges as well as the liquids and secretions suctioned during medical procedures. The disposable cartridges are made of a material that does not emit harmful substances during combustion. The combustion or other utilization procedures should be carried out in accordance with the internal regulations provisions of a particular health service unit.

Catalogue №	A-OG1000	A-OG1000Z	A-OG2000	A-OG2000Z
Name	<b>OG1000</b>	<b>OG1000Z</b>	<b>OG2000</b>	<b>OG2000Z</b>
Capacity	1l	1l	2l	2l
Gelling agent	no	yes	no	yes